




ANAEROBIC GLOVE BOX



01 ACCURACY

- Achieve anaerobic conditions in 20-30 minutes
- Can maintain a stable anaerobic environment (< 5 PPM) for weeks
- Save gas after achieving an anaerobic environment
- Precise Temperature Control
 - ✓ Range of Temperature Control: 5°C above ambient up to 45°C, in 0.1°C increments
- Advanced Airway Design
 - ✓ Ensures a high level of homogeneity in the workstation's gas distribution



- Easy Access Operation System

✓ The sleeves are made of a soft, skin-friendly material and can access the workstation easily and quickly without vacuum



Easy-access Operation System



Direct-access Operation System

- Direct-access operation system: Optional sleeve vacuum system
- Interlock can be quickly purged to $< 0.1\% \text{ O}_2$ within 75 seconds
- USB data logging system holds 3 months of continuous data storage

02 CUSTOMISED PRODUCTS



Anaerobic Glove Box with a Larger Workspace YY-TME



Anaerobic Glove Box with a Larger Workspace YY-TME Plus

Tel: 8610-88693537

Email: sales@maworde.com

www.maworde-biotech.com

