


HYPOXIA WORKSTATION FOR FISH RESEARCH



Specification

- A large fish tank with a mixer is integrated. The fish tank is connected to an external water pipe with a switch, which make it easy to change water.
- 8 extra small tanks with water pumps connected to the large fish tank
- Sleeve ports hand-direct access operating system facilitates, Avoiding Function cross-contamination and skin allergies caused by rubber gloves.
- Standard built-in interlock allows ease of transferring samples into the workstation
- O₂ control range:
0.1-23.0%, in increments of 0.1%.
- USB data logging system:
Set per minute, with time, date, and O₂ (set/actual)
- Optional temperature control, dissolved oxygen detector, etc.
- Providing a customised size based on your needs

Tel: 8610-88693537

Email: sales@maworde.com

www.maworde-biotech.com

