

# TRI GAS MIXER (QTK-A)

Mix 100% N<sub>2</sub>, CO<sub>2</sub> and O<sub>2</sub> to get various concentrations of O<sub>2</sub> and CO<sub>2</sub>.



## Features & Benefits

- **Precise Gas Control:**

- O<sub>2</sub> control range: 0.1% to 23.0%, in 0.1% increments
  - CO<sub>2</sub> control range: 0.1% to 30.0%, in 0.1% increments
  - Programmable four-stage hypoxic cycle

- **Data Logging System:**

- USB data logging system holds 3 months of continuous data storage: time, date, O<sub>2</sub> (set/actual), CO<sub>2</sub> (set/actual).



Provide fixed proportion of N<sub>2</sub>, CO<sub>2</sub> and O<sub>2</sub> for imaging system.



Provide fixed proportion of N<sub>2</sub>, CO<sub>2</sub> and O<sub>2</sub> for Microplate Reader or high content imaging system.

Tel: 8610-88693537

Email: [sales@maworde.com](mailto:sales@maworde.com)

[www.maworde-biotech.com](http://www.maworde-biotech.com)

